


Principle of Confocal Microscopy and its applications


Confocal Microscopy



To overcome some limitations of traditional wide-field fluorescence microscopes.

Confocal Microscopy is now frequently known as

- **confocal laser scanning microscopy (CLSM)**
 - **laser confocal scanning microscopy (LCSM)**,
- is an optical imaging technique for **increasing optical resolution and contrast** of a micrograph by means of using a **spatial pinhole to block out-of-focus light** in image formation. (From wiki)



Concept of Confocal Microscopy

Confocal Microscopy-Spatial Resolution

For a Confocal Microscope, the pinhole in front of sensor should be smaller than the diameter of Airy diffraction image formed by the object lens. The image of Airy diffraction pattern would be trimmed to region around the center.

Resolution of a optical microscope

$$\text{Lateral} = \begin{cases} \frac{\lambda}{2 \times NA}, \text{Abbe Equ.} \\ \frac{0.61\lambda}{NA}, \text{Rayleigh limit} \\ \frac{1.22\lambda}{NA_{obj} + NA_{cond}}, \end{cases}$$

$$\text{Depth} = \frac{2\lambda}{NA^2}, \text{Abbe}$$






Figure 1

Resolution of a optical microscope


Lateral Resolution = $0.37\lambda/NA$

$$\text{Depth Resolution} = \begin{cases} \sqrt{\left(\frac{0.88\lambda}{1 - \sqrt{1 - NA^2}}\right)^2 + \left(\frac{\sqrt{2}d}{NA}\right)^2}, \text{Geometrical-optical confocality} \\ \frac{0.64\lambda}{1 - \sqrt{1 - NA^2}}, \text{Wave-optical confocality} \end{cases}$$


DOI 10.1007/978-3-642-12012-1




Types of Confocal Microscopes



Laser Scanning Confocal Microscope










Disc Scanning Confocal Microscope





Programmable Array Scanning Confocal Microscope


Laser Scanning Confocal Microscope- Scanning Mechanisms





Confocal Fluorescence Microscopy



Confocal Microscopy:
higher z-resolution and reduced
out-of-focus-blur



Optical Section



Spinning disk confocal


Pros:

- Accelerating the scanning process
- Fast – multiple points are illuminated at once
- Photon efficient – high QE of CCD
- Gentler on live samples – usually lower laser power

Cons:

- multiple fluorophore colors might not be generated because of PMT is replaced by CCD sensor
- Cross talk problem would degrade depth resolution
- Fixed pinhole – except in swept-field
- Small field of view (usually)

CCD and the pinhole disk must be synchronized



Yokogawa Spinning Disk Scanning Unit









Figure 1










Programmable Array Scanning Confocal Microscope



SCANNING VOL. 38, 234–239 (2016)




DMD(Digital micromirror device)



- Switching Speed of the ferroelectric liquid crystal is $40 \mu s$
- FLCos device
 - CMOS substrate/ individual pixels as mirror
 - Cover front glass coted with a transparent electrode
 - middle layer is the ferroelectric liquid crystal
- On/Off by switching polarization angle:
zero angle polarization state is manufactured as the black point of the device, while the 90° polarization state is the white point

Chromatic (dispersion) confocal microscopy (CDCM)




Lateral Scanning Mechanism

- Laser Scanning Confocal Microscope
- Disc Scanning Confocal Microscope
- Programmable Array Scanning Confocal Microscope

Vertical Scanning—Optical Section



- Moving the Objective lens



Vertical Scanning is replaced by Chromatic Aberration → the optical head entirely static, without any spurious vibration generated by an internal mechanism




Recalling Chromatic Aberration, different wavelength would focus at different focal point along the optical axis

Line-type chromatic confocal sensors







https://www.focalspec.com/wp-content/files/LCI_animation_14112018.mp4

Line-type chromatic confocal sensors




The tilted focal plane is formed from the tilted pinhole mask by an objective lens. To create a parallel depth scanning scheme, the specimen is scanned through the stationary stack of diffraction foci in a stepwise fashion or in a continuous approach

Area Type Chromatic Confocal Microscope



One Shot




Area Method

Chromatic confocal detection for high speed micro-topography measurements

Chromatic confocal matrix sensor

Separate illumination and detection pinholes






<http://dx.doi.org/10.1364/AO.53.007634>

<http://dx.doi.org/10.1364/AO.54.004927>


Chromatic confocal matrix sensor

Single pinhole array used in double pass









Dual Sensors





Dual Sensors



Determined by Ratio of Light Intensity Detected at PMT1 and PMT2
In front of PMT1, there is an optical filter placed in front of the sensor



Dual Sensors



doi:10.1088/0957-0233/25/12/125403

